UNIVERSITY OF MIAMI



Rosenstiel School of Marine, Atmospheric, and Earth Science Tritium Laboratory 4600 Rickenbacker Causeway Miami, FL 33149-1031 P: 305-421-4100 F: 305-421-4112 tritium@miami.edu

Tritium Laboratory 23 February 2023

SWAB REPORT # 1053

SWAB DATE: 17 February 2023

*R/V Atlantic Explorer*UNOLS Radioisotope Van #625.5.02

James D. Happell

Distribution: SWAB Committee Capt. Rick Verlini Rod Johnson Angie Knapp Tim Deering

COMMENTS TO SWAB REPORTS

The LSC is now a Quantulus GCT 6220, with the SWAB counting assay having background cpm of 0.3 & 1.2 for ³H & ¹⁴C. This replaces an LSC with background cpm of 1.6 & 5.5 for ³H & ¹⁴C.

All samples are counted for 60 minutes, the instrument background is subtracted, and activities are reported in dpm/m². Bucket blank activities are not subtracted. Counting errors (2 standard deviations) are also reported in dpm/m². An error larger than the activity indicates that the activity is not significantly different from zero. All activities significantly above background will be in **bold**.

Criteria for SWAB Results

Category	3 H (dpm/m 2)	14 C (dpm m ²)	Recommendations
A	< 500	<50	No action
B*	500-10,000	50-10,000	Needs cleaning before any natural tracer work. Decks in radiation vans with activities
			above 1000 dpm/m ² should be cleaned.
C**	10,000-100,000	10,000-50,000	Must be cleaned before any use.
D***	>100,000	>50,000	May be a health hazard. Notify local radiation safety official.

Note: ¹⁴C and ³⁵S have peak energies of 156 and 167 KeV, respectively; thus ³⁵S will be registered as ¹⁴C by our counting techniques. Categories A, B and C are not a health hazard.

<u>Recommended Cleaning Proceedure</u> Wearing ordinary household rubber gloves:

Disposal of Cleaning Materials (gloves, sponges, etc)

Categories A & B dispose as ordinary garbage, C & D contact your institution's radiation safety office.

Note: If category C or D is encountered, we try to notify the insitution promptly by phone or email.

³H: Wash and scrub with radioactive cleanup detergent such as COUNT-OFF (50 ml COUNT-OFF to 4 liters of water), using sponges to distribute solution and reabsorb it.

¹⁴C: Wash with 1% sulfuric or 2% hydrochloric (muriatic) acid with good ventilation (will dissolve carbonates, releasing ¹⁴CO₂). Follow up with wash as if for ³H.

REPORT FOR SWAB # 1053

LOCATION: St. Petersburg, FL VESSEL: *R/V Atlantic Explorer*


DATE: 16 February 2023 TECHNICIAN: Jim Happell

Sample # Sample Identification	³ H dpı	³ H dpm/m ²		¹⁴ C dpm/m ²		
	activity	error	activity	-	error	
1 1st Vial Bkgnd	0	± 0	0	±	0	
2 Initial bucket blank	-12	± 50	-12	±	40	
Forward Lab (Figure 1)						
3 Forward benchtop	-14	± 56	-5	\pm	17	
4 Port benchtop forward of sink	-14	± 59	-21	\pm	74	
5 Deck inside starboard entrance	-16	± 65	-30	\pm	104	
6 Port sink area	2	\pm 13	-9	\pm	33	
7 Center benchtop	-19	± 79	0	\pm	2	
8 Port benchtop aft of sink	-13	± 53	-10	\pm	34	
9 Deck outside Enviro Room	7	± 40	-21	\pm	73	
10 Deck in front of sink	-13	± 51	-7	±	26	
Main Lab (Figure 1)						
11 Deck inside forward entrance	1	± 6	-13	\pm	47	
12 Port sink area	-29	± 86	-24	\pm	83	
13 Benchtop inside laminar flow hood	-26	± 77	-8	\pm	29	
14 Deck inside aft entrances	-8	\pm 0	-10	\pm	34	
15 Inside Whirlpool freezer	23	\pm 67	-28	\pm	96	
16 Inside Frigidaire freezer	-24	\pm 71	-15	\pm	51	
17 Center benchtop	-40	\pm 118	0	\pm	37	
18 Port benchtop aft of sink	-19	± 75	-13	\pm	46	
19 Port benchtop forward of sink	-7	\pm 39	-14	\pm	50	
20 Aft startboard benchtop	-12	± 48	-21	±	72	
Main Lab (Figure 1)						
21 Inside fume hood	9	\pm 51	-26	\pm	90	
22 Forward benchtop	-27	\pm 80	-17	\pm	58	
23 Inside -80°C freezer #2	-43	± 126	-16	\pm	57	
24 Port benchtop forward of sink	-42	± 125	-2	\pm	19	
25 Port sink area	-21	± 62	2	\pm	23	
26 Deck inside starboard aft entrance	2	± 11	-15	\pm	53	
27 Forward section of center benchtop	-22	± 65	-4	\pm	15	
28 Inside DEAD Cospolich refrigerator	-21	± 63	-19	±	66	
29 Inside LIVE Cospolich refrigerator	-28		-17	\pm	59	

Sample # Sample Identification	³ H dpm/m ²	¹⁴ C dpm/m ²		
	activity error	· activity error		
30 Deck inside forward entrance	-14 ± 59	-2 ± 15		
31 Deck in front of port sink	10 ± 42	-7 ± 26		
32 Inside -80°C freezer #1	-19 ± 76	-21 ± 72		
33 Inside DEAD Cospolich freezer	-24 ± 70	1 ± 7		
34 Inside LIVE Cospolich freezer	-22 ± 66	5 -5 ± 17		
Radioisotope Van (Figure 2)				
35 Deck outside van entrance	-36 ± 107	6 ± 19		
36 Sink area	-2 ± 13	-3 ± 10		
37 Benchtop adjacent to sink	99 ± 36	11 ± 10		
38 Benchtop adjacent to fume hood	6 ± 38	-4 ± 15		
39 Inside refirgerator	80 ± 30	39 ± 14		
40 Benchtop adjacent to LSC	-10 ± 42	$2 4 \pm 15$		
41 Benchtop across from sink	26 ± 32	-5 ± 18		
42 Center Deck	280 ± 56	48 ± 13		
43 Deck inside entrance	81 ± 38	29 ± 13		
44 Final bucket blank	-44 ± 129	-1 ± 9		

Comments

Please note that the error reported for each isotope is the two-standard deviation counting error. Reports may now contain values less than zero. Decay counting background samples will be distributed about the background vial, which means that negative values are possible. In the past we rounded the negative values to zero. Values are only significantly above background when they are positive and larger than the error. Please note that we are now using a Quantulus 6220 LSC which counts very near natural background. While the cleanup standards have not changed all values above background will now be in bold. There was no isotope contamination that requires cleaning detected in samples from the ship or from the rad van.






Figure 2 SWAB #1053 17 February 2023