



Rosenstiel School of Marine, Atmospheric, and Earth Science Tritium Laboratory 4600 Rickenbacker Causeway Miami, FL 33149-1031 P: 305-421-4100 F: 305-421-4112 tritium@miami.edu

Tritium Laboratory 27 August 2025

SWAB REPORT # 1126

SWAB DATE: N/A

University of Nebraska

Dr. James D. Happell Associate Research Professor

Distribution: SWAB Committee David Manning The LSC is now a Quantulus GCT 6220, with the SWAB counting assay having background cpm of 0.3 & 1.2 for ³H & ¹⁴C. This replaces an LSC with background cpm of 1.6 & 5.5 for ³H & ¹⁴C.

All samples are counted for 60 minutes, the instrument background is subtracted, and activities are reported in dpm/m². Bucket blank activities are not subtracted. Counting errors (2 standard deviations) are also reported in dpm/m². An error larger than the activity indicates that the activity is not significantly different from zero. All activities significantly above background will be in **bold**.

Criteria for SWAB Results

Category	3 H (dpm/m 2)	14 C (dpm m 2)	Recommendations		
A B*	<500	<50	No action		
Β.	500-10,000	50-10,000	Needs cleaning before any natural tracer work. Decks in radiation vans with activities		
			above 1000 dpm/m ² should be cleaned.		
C**	10,000-100,000	10,000-50,000	Must be cleaned before any use.		
D***	>100,000	>50,000	May be a health hazard. Notify local radiation safety official.		

Note: ¹⁴C and ³⁵S have peak energies of 156 and 167 KeV, respectively; thus ³⁵S will be registered as ¹⁴C by our counting techniques. Categories A, B and C are not a health hazard.

<u>Recommended Cleaning Proceedure</u> Wearing ordinary household rubber gloves:

Disposal of Cleaning Materials (gloves, sponges, etc)

Categories A & B dispose as ordinary garbage, C & D contact your institution's radiation safety office.

Note: If category C or D is encountered, we try to notify the insitution promptly by phone or email.

³H: Wash and scrub with radioactive cleanup detergent such as COUNT-OFF (50 ml COUNT-OFF to 4 liters of water), using sponges to distribute solution and reabsorb it.

¹⁴C: Wash with 1% sulfuric or 2% hydrochloric (muriatic) acid with good ventilation (will dissolve carbonates, releasing ¹⁴CO₂). Follow up with wash as if for ³H.

REPORT FOR SWAB # 1126

LOCATION: Omaha, NE DATE: N/A

VESSEL: particulates on filter paper TECHNICIAN: David Manning

Sample # Sample Identification	³ H dpm			¹⁴ C dpm		
	activity		error	activity		error
1 1st Vial Background	7.4	土	1.2	-2.3	土	1.6
2 1 Filter paper with stream particulates	11.6	\pm	1.7	-2.2	\pm	1.6
3 2 Filter paper with stream particulates	11.6	\pm	1.7	-2.0	\pm	1.4
4 3 Filter paper with stream particulates	10.9	\pm	2.3	-0.6	\pm	0.4
5 4 Filter paper with stream particulates	15.3	\pm	2.7	-0.2	\pm	0.1
6 5 Filter paper with stream particulates	9.8	\pm	1.7	-1.6	\pm	1.1
7 6 Filter paper with stream particulates	10.6	\pm	1.6	-1.9	\pm	1.4
8 7 Filter paper with stream particulates	11.0	\pm	2.0	-0.8	\pm	0.5
9 8 Filter paper with stream particulates	13.5	\pm	2.7	0.1	\pm	0.0
10 9 Filter paper with stream particulates	11.9	\pm	2.4	-0.4	\pm	0.3
11 10 Filter paper with stream particulates	12.4	\pm	1.8	-1.7	\pm	1.2
12 11 Filter paper with stream particulates	13.9	\pm	2.5	-0.1	\pm	0.0
13 12 Filter paper with stream particulates	13.0	\pm	1.9	-2.0	\pm	1.4
14 13 Filter paper with stream particulates	11.1	\pm	2.0	-0.7	\pm	0.5
15 14 Filter paper with stream particulates	15.6	\pm	2.9	0.0	\pm	0.0
16 15 Filter paper with stream particulates	14.0	\pm	2.5	0.1	\pm	0.0
17 16 Filter paper with stream particulates	12.1	土	2.4	0.1	\pm	0.0
18 17 Filter paper with stream particulates	12.9	\pm	2.4	-0.4	\pm	0.3
19 18 Filter paper with stream particulates	11.8	±	1.9	-1.2	±	0.8

Comments

Please note that the error reported for each isotope is the two-standard deviation counting error. Reports may now contain values less than zero. Decay counting background samples will be distributed about the background vial, which means that negative values are possible. In the past we rounded the negative values to zero. Values are only significantly above background when they are positive and larger than the error. Please note that we are now using a Quantulus 6220 LSC which counts very near natural background. The clean-up standards do not apply in the case. Since we are counting filter paper with particulates and not the typical water and countoff solution that has been spread on a horizontal surface, we are reporting the results in dpm rather than dpm/m². There is no indication on the LSC of ¹⁴C contamination